



# TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

## Inorganic Carbon Standard Worksheet

Date, Name of preparer		line 1
Lot of Na <sub>2</sub> CO <sub>3</sub> used in prep.		line 2
Weight of Na <sub>2</sub> CO <sub>3</sub> + weigh boat		line 3
Weight of weigh boat		line 4
Weight of Na <sub>2</sub> CO <sub>3</sub>		line 5

Table 1. **Preparation of the Stock Na<sub>2</sub>CO<sub>3</sub> standard**

Concentration of stock solution		line 6
	nominal 500 ppm	

Table 2. **Preparation of the Dilute Na<sub>2</sub>CO<sub>3</sub> standard**

Concentration of dilute solution	nominal 25.0 ppm	line 7
----------------------------------	------------------	--------



# TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. [MICLAB 115](#))

## Organic Carbon Standard Worksheet

Date, Name of preparer		line 1
Lot of KHP used in prep.		line 2
Weight of KHP + weigh boat		line 3
Weight of weigh boat		line 4
Weight of KHP		line 5

Table 1. **Preparation of the Stock KHP standard**

Concentration of stock solution		line 6
	nominal 2000 ppm	

Table 2. **Preparation of the Dilute KHP standard**

Concentration of dilute solution	nominal 25.0 ppm	line 7
----------------------------------	------------------	--------



# TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

## Calibration Worksheet for Autosampler Calibration

Date, Name of Analyst \_\_\_\_\_

### Measurements of Standards

Average TC of 3rd DI blank		=TC <sub>DI blank</sub>
Average TC of 3rd TC(KHP) standard		=TC <sub>ave,TC</sub>
Average TC of 3rd IC (Na <sub>2</sub> CO <sub>3</sub> ) standard		=TC <sub>ave,IC</sub>
Average IC of 3rd IC (Na <sub>2</sub> CO <sub>3</sub> ) standard		= IC <sub>ave,IC</sub>
TC value of TC(KHP) standard from formulation		= TC <sub>std</sub>
Adjusted TC value of the TC(KHP) standard = TC <sub>std</sub> + TC <sub>DI blank</sub>		= TC <sub>adj,TC</sub>

### Calculation of Calibration Constants

Current TC calibration constant		=TC <sub>cal,old</sub>
---------------------------------	--	------------------------

New TC calibration constant	$= \frac{TC_{cal,old} \cdot TC_{adj,TC}}{TC_{ave,TC}}$	=TC <sub>cal,new</sub>
-----------------------------	--	------------------------

Current IC calibration constant		=IC <sub>cal,old</sub>
---------------------------------	--	------------------------

New IC calibration constant	$= IC_{cal,old} \cdot \frac{TC_{adj,TC}}{TC_{ave,TC}} \cdot \frac{TC_{ave,IC}}{IC_{ave,IC}}$	=IC <sub>cal,new</sub>
-----------------------------	--	------------------------

### Verification of Calibration

Average TC of 3rd TC (KHP) standard		=TC <sub>ave,chk,TC</sub>
Average TC of 3rd IC(Na <sub>2</sub> CO <sub>3</sub> ) standard		=TC <sub>ave,chk,IC</sub>
Average IC of 3rd IC (Na <sub>2</sub> CO <sub>3</sub> ) standard		=IC <sub>ave,chk,IC</sub>
New TC <sub>adj,TC</sub> = TC(KHP) <sub>std</sub> + average TC 3 <sup>rd</sup> DI Blank =		
TC Std % difference = $\frac{(New\ TC_{adj,TC} - TC_{ave,chk,TC}) \times 100}{New\ TC_{adj,TC}}$ =		(NMT 3%)
IC Std % difference = $\frac{(TC_{ave,chk,IC} - IC_{ave,chk,IC}) \times 100}{TC_{ave,chk,IC}}$ =		(NMT 3%)



# TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

## 500 ppb Sucrose Standard Worksheet

Date, Name of Preparer	line 1
Lot No. of Sucrose used	line 2
Carbon Assay (% value / 100)	line 3 nominal 0.421
Gross Weight (Sucrose + Boat)	line 4
Tare Weight of Weigh Boat	line 5
Nett Weight of Sucrose	line 6
Carbon Conc. of Stock Solution (ppb C = $\frac{\text{g [line 6]} * \text{ca [line 3]} * 10^6}{1 \text{ L}}$ )	line 7 nominal 25,000 ppb C
Carbon Conc. of Working Solution ([line 7] ÷ 50)	line 8 nominal 500 ppb C



# TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. [MICLAB 115](#))

## 500 ppb Benzoquinone Standard Worksheet

Date, Name of Preparer	line 1
Lot No. of Benzoquinone used	line 2
Carbon Assay (% value / 100) nominal 0.666	line 3
Gross Weight (Benzoquinone + Boat)	line 4
Tare Weight of Weigh Boat	line 5
Nett Weight of Benzoquinone nominal 0.075g	line 6
Carbon Conc. of Stock Solution (ppb C = $\frac{\text{g [line 6]} * \text{ca [line 3]} * 10^6}{1 \text{ L}}$ ) nominal 50,000 ppb C	line 7
Carbon Conc. of Working Solution ( $[\text{line 7}] \div 100$ ) nominal 500 ppb C	line 8



# TOC ANALYSER CALIBRATION WORKSHEETS

(Ref. MICLAB 115)

## Suitability Verification Worksheet

Date and Name of Analyst		line 1
--------------------------	--	--------

Calculated Concentration of Sucrose Standard		line 2
--	--	--------

Calculated Concentration of Suitability Standard		line 3
--	--	--------

Average TOC of Water Blank		line 4
----------------------------	--	--------

Average TOC of Sucrose Standard		line 5
---------------------------------	--	--------

Average TOC of Suitability Standard		line 6
-------------------------------------	--	--------

<b>Response Efficiency</b> Response Efficiency (%) = $\frac{(TOC_{suitability} - TOC_{blank})}{(TOC_{sucrose} - TOC_{blank})} * \frac{C_{sucrose}}{C_{suitability}} * 100$	between 85% and 115%	line 7
---	----------------------	--------