




Labelling of Settle (Fallout) Plates

(Ref. [MICLAB 095](#))

Note: Before logging settle plates into the sampling cabinet always check the following –

1. The plate is an **NA Agar** plate and is within **expiry**.
2. The plates have no growth on them.
3. **ALL** of the following information is **neatly** written on the plate -

- 
- BPN
 - Material Code Number
 - Date of Filling
 - Machine Name and Number
 - Start time of Exposure
 - Finish time of Exposure
 - Position on machine i.e. Left, Right side

Failure to include any of the above mentioned information will result in the creation of a **DR**.